

Carestream - CBCT - saving screenshots into TDO

342 [TDO KB](#) November 21, 2025 [CBCT](#), [Troubleshooting](#) 28.04K

Overview


The article below reviews the steps to save a Carestream CBCT screenshot into TDO. If you are experiencing issues with this function, it is also advisable to review the following steps to ensure everything is set up correctly. Once a CBCT volume is acquired, users can take as many screenshots as necessary. If the volumes are open from TDO, screenshots taken on Kodak 3D Viewer will automatically be saved in TDO.

CS 3D Imaging Software version 3.10.38 or newer is recommended

Detailed Instructions

The following steps must be manually applied to each computer that is intended to be used to capture Carestream CBCT screenshots and saved into TDO.

- Within the Carestream Viewing Application, select **Settings** (shaped like a gear)
- From the **Preferences** options, select **Screenshots** (shaped like a camera)



- If **Save a copy in the specified folder** is checked off, uncheck it
- Underneath the **Screenshot format** section, select **TIFF**
- Click **Apply** to save the changes
- Close the Carestream application

Saving Carestream screenshots into TDO

There are two (2) methods of taking screenshots that can be saved automatically in TDO:

Method 1: immediately after CBCT acquisition

- From the TDO **Acquire** window, select the **CBCT** button to take a new CBCT
- After the CBCT has been acquired, within the Carestream application, set up what you'd like to screenshot and click **Screenshot** (camera icon)
- The screenshot will automatically save into the TDO **Image Organizer**

Method 2: existing CBCT

- From the TDO **Image Organizer**, open the CBCT volume that you'd like to screenshot, either **double-clicking** or **right-clicking** and selecting **Open Volume**
- From the Carestream application, set up what you'd like to screenshot and click **Screenshot** (camera icon)
- The screenshot will automatically save into the TDO **Image Organizer**

Online URL: <https://kb.tdo4endo.com/article.php?id=342>