

J Morita (i-Dixel) CBCT Integration

715 [TDO KB](#) November 20, 2025 [CBCT](#) 21.07K

Overview

TDO has integrated J Morita's i-Dixel CBCT units for a direct connection to their 3D Imaging software. Patients' files are linked, allowing easy access between programs in one window, which eliminates file duplication. TDO users can now seamlessly integrate a Morita 3D unit into their practice with full compatibility.

Requirements

i-Dixel Version 1.817 or later (JMPatientinfo.dll version 1.0.0.8 or later). If on an older version, please contact J Morita technical support to update the i-Dixel version.

We do not currently support iDixel Web.

TDO Version 11.111 or later

Note: Your IT must ensure that all Windows users have full read/write permission on the following directories:

Registry folders:

HKEY_LOCAL_MACHINE\SOFTWARE\JMorita

HKEY_LOCAL_MACHINE\SOFTWARE\WOW6432Node\JMorita

HKEY_LOCAL_MACHINE\SOFTWARE\WOW6432Node\ODBC

C drive:

C:\3dx_data

C:\3dx_settings

C:\JmoritaMFG

C:\Program Files (x86)\3dx

C:\Program Files (x86)\3dxcom

C:\Program Files (x86)\3dxdbmnt

C:\Program Files (x86)\JMorita

Configuration Version 11.1-12.2

1. Open Administrative > Edit Control Table.
2. All computers must have a unique name. If you have two or more computer with the same name they must be renamed. If the duplicate computer name is no longer in use please contact support

to remove from your database.

3. Select the Workstation Settings tab and fill in the CBCT Settings section.

The Default CBCT Image Category will be used to save the Morita Volumes and screenshots under the specified category. You can create a new Image Category just for CBCT's by going to Setup > Program Setup > Image Categories, it is typically Morita CBCT.

CBCT Software: Morita

Morita Server: Varies, but this is usually 3DX-SERVER. If you're not sure, open i-Dixel preferences, and go to the Database section.

Morita Database: Varies, but is usually 'master' or 'user_net'

Morita UID: Varies, but is usually '1000' or '500'

Morita PWD: 1

Morita Screenshot Folder: C:\JMoritaMFG\3dx\temp\TDO (make sure that this folder exists and that JMC Tech Support has configured 3rd party screenshots to go to this folder)

Control Table Information

Practice Information **Application Settings** **Network Settings** **Workstation Settings**

Digital Radiography
Radiography Software: [Dropdown]
Radiography Path/Source: [Text Box]

Remote Camera
Software: [Breeze Systems]
Remote Camera Software Path: [c:\Program Files\BreezeSys\DSLR Remote Pro]
☐ Autosave Camera Image ☐ Autoload Remote Software

CBCT Settings
CBCT Software: [Morita] Default CBCT Image Category: [CBCT]
Morita Server: [3dx-server] Morita Database: [master] Morita UID: [1000] Morita PWD: [1]
Morita Screenshot Folder: [C:\JMoritaMFG\3dx\temp\TDO]

Image Location
Flash Card Reader Path: [Text Box]

Patient Privacy
☐ Open Application Showing Calendar First

Tooth Models
Tooth Model Path: [C:\Program Files\Brown & Herbranson\3D Tooth Atlas]

Digital Signature
Signature Type: [Leadtools (Default)]

Notes Encryption
☐ Allow Notes Encryption

Inactivity Duration Before Logoff
Minutes: [20]


Pin Pad Settings
Pin Pad Type: [Dropdown] Connection Method: [Dropdown] Com Port: [Text Box]
☐ Canadian EMV

Close


Version 12.311b+

Note: For CBCT Unit configuration in 12.311b or newer, please click [here](#).

1. Open Setup > Site Configuration > Device Configuration
2. Select the Workstation you need to configure
3. Select Device: Morita, then click New
4. Check the '**CBCT is attached to this computer!**' box; **CBCT Default Category:** CBCT (or whatever your practice uses); **Server:** Usually 3DX-SERVER or 3DX-CAPTURE; **Database:** usually user_net or master; **Username:** 1000; **Password:** 1; **Screenshot:** C:\JMoritaMFG\3dx\temp\TDO (make sure that this folder exists and that **JMC Tech Support has configured 3rd party screenshots** to go to this folder)



5. Click the Save button (floppy disk icon, upper right-hand corner)




Acquiring a 3D Volume

TDO will connect to the i-Dixel application and automatically check if the patient already exists or not. This information is checked by comparing the TDO Patient ID and i-Dixel patient ID, E.g.: TDO Patient ID 234 = i-Dixel Patient ID: TD00000234. If the Patient does not exist in the i-Dixel database TDO will automatically create it. The i-Dixel application requires the following information to create a patient, make sure the following patient information is entered in TDO:

1. First Name
2. Last Name
3. Gender
4. Birth Date

Steps to Acquire

1. On the Main Page of TDO, click on Acquire
2. Click the CBCT button, this will launch i-Dixel



3. Acquire the CBCT as you normally would in i-Dixel

Online URL: <https://kb.tdo4endo.com/article.php?id=715>