

Entering Anesthesia

113 [TDO KB](#) November 26, 2025 [New User resources](#) 4.80K

Overview

This article outlines how to record anesthesia that the patient was given in TDO.

Detailed Instructions

There are buttons to open the Anesthesia page on both the **Treatment** page and the **Prescription** page.

The screenshot displays the TDO software interface for entering anesthesia. At the top, there is a 'History' tab with columns for Date, Dosage, Frequency, Medication Prescribed, Reason for Taking, Stopped Taking, Date Last Taken, and Anesthesia. Below this is a large blue rectangular area for data entry. At the bottom, there is a form with various fields and buttons. The 'Anesthesia' button is highlighted in yellow. Other buttons include 'Complete Group', 'Occlusion Reduction', 'Fractures', and 'XRay Data'. There are also checkboxes for 'Email Interim Report' and 'Fractures'. A 'Cases filter' section on the right has radio buttons for 'All', 'Group', 'Tooth', and 'Pick'.

The dental assistant fills out the anesthesia information as the doctor administers the anesthetic. Recorded are the kind of anesthesia, the type of injection used, the amount given, and the date and time administered.

You can customize all of these lists by pressing the + button or by right-clicking and selecting **Edit List**.

There is also a way to auto-populate this form if you routinely give the same kind and amount of anesthesia. The Anesthesia Autofill feature fills in the anesthesia, type, quantity, and date all at once. You can customize this auto-fill routine from the [Anesthesia Autofill options](#). To use an autofill, select it from the combo box and click the arrow button.

Anesthesia Autofills Options

Anesthesia Autofills

Name: Carbo
Anes: 2% Lido no epi
Type: MB
Qty: 1 Carp

Name: Lido epi
Anes: 2% Lido 1:50,000 epi
Type: I
Qty: 1 Carp

Name: Marc-I
Anes: Marcaine 1:200,000
Type: I
Qty: 1 Carp

Name: Marc-blk
Anes: Marcaine
Type: MB
Qty: 1 Carp

Record: 1 of 7 | No Filter | Search

There is also a bulk anesthesia entry if you frequently give the same multiple injections - for example, a mandibular block and a long buccal injection. These are customized in Anesthesia Explosion Groups. To use an explosion autofill, select it from the combo box (the combo box to the right of the individual autofill) and click the arrow button.

Anesthesia Explosion Groups

Name: test

Anesthesia entries to be added when this explosion code is selected:

▶ 1 Carp	MB	2% Lido 1:100,000 epi
1/2 Carp	STBDT	2% Lido 1:100,000 epi
*		

Record: 1 of 2 | No Filter | Search

The other part of this form documents the Nitrous Oxide parameters, if needed. The form is a pop-up and is self-explanatory.

Online URL: <https://kb.tdo4endo.com/article.php?id=113>