

Adding a watermark to images in TDO

691 [TDO KB](#) February 26, 2026 [Organizer](#), [Organizer](#), [Site Configuration](#) 6.56K


Overview

This article will show you how to add a custom watermark to images or radiographs in TDO using a .PNG version of your logo. If you do not currently have a .PNG version of your logo, reach out to the graphic designer who created your logo.


Note: This feature requires TDO Version 12.205

Detailed Instructions

1. Go to the **Setup > Site Configuration** and select **Edit Control Table**.
2. Click on the **Practice Information** tab, and in the lower left corner, you will see the **Logo for Reports** section (see image)
3. Click the **Browse...** button directly below that, and select the .PNG version of your logo.



4. Once the image is loaded, it will show beneath the **Browse...** button.



5. Make sure that the position of the watermark is where you would like it to be.

Online URL: <https://kb.tdo4endo.com/article.php?id=691>

