

Maximizing the TDO Window in TDO 12.338+

764 [TDO KB](#) April 22, 2025 [Site Configuration](#) 6.02K

Overview

In TDO 12.338 and higher, it is possible to enable a setting that will allow you to maximize the main window so that it occupies the whole screen.

Detailed Steps

- Open TDO, and log in.
- Go to Setup > Site Configuration > Edit Control Table
- Select the Workstation Settings tab.
- Check the box that says Allow TDO to be Full Screen

The screenshot shows the 'Control Table Information' dialog box with the 'Workstation Settings' tab selected. The 'Windows Placement' section at the bottom right contains a checkbox labeled 'Allow TDO to be Full Screen' which is checked. A red checkmark is drawn over this checkbox. Other visible settings include 'Digital Radiography' (Carestream 6000 Series), 'Scan Twain Settings' (Default selected), 'Remote Camera' (Breeze Systems), 'Image Location' (Flash Card Reader Path), 'Patient Privacy' (Open Application Showing Calendar First), 'Tooth Models' (C:\Program Files\Brown & Herbranson\3D Tooth Atlas), 'Digital Signature' (Leadtools (Default)), 'Inactivity Duration Before Logoff' (999 minutes), and 'Credit Card Settings' (Unknown company).

- Close the Control Table. You will now be able to click the Maximize button (between Minimize and Close) to make the main TDO window full-screen.

Online URL: <https://kb.tdo4endo.com/article.php?id=764>