

Carestream - Trophy - integration

860 [TDOKB1](#) November 18, 2025 [Radiography](#), [Site Configuration](#) 11.49K

Overview

This is for setting up the trophy integration to use RVG with Carestream Imaging software.

Detailed Instructions

This is the best adjustment for offices that have issues with Carestream's SDK.

Setup > Site Configuration > Edit Control Table > Workstation Settings

The screenshot shows the 'The Digital Office' software interface. The 'Setup' menu is open, and 'Site Configuration' is selected, leading to the 'Edit Control Table' dialog. The patient information form is visible, showing details for Michael McNair. The form includes fields for personal information, contact details, and medical history. A 'Remove Patient' button is located at the bottom left of the form.

Site Configuration Menu:

- Edit Control Table
- Edit Doctor Information
- Electronic Claims Setup
- Add/Edit User Names and Passwords
- Add/Edit Practice and Doctor Emails
- Employee Information
- Equipment Information
- Software Information
- Device Configuration

Patient Information Form:

Personal Information:

- Last Name: McNair
- First Name: Michael
- Phone: 202252
- Salutation: Mr.
- Address1: 7732 Middle Valley Drive
- Address2:
- City: Springfield
- State: VA
- Zip: 22153
- Phone(Home): (703) 644-6326 x
- Phone(Work): (202) 252-6494 x
- Phone(Mobile): (703) 850-4659 x
- Phone(Fax):
- SS#: *****
- Birthdate: 7/3/1951
- Age: 69
- Gender: Male

Medical History:

- Referred Date: 6/1/2011
- ID: 4
- Chart #:

Referred by:

- Referred Practice:
- Gen Dentist:
- Additional Doctor:
- Occupation:
- Employer Name: BAE Systems
- Employer Address:
- Emerg. Contact: Glennis McNair
- Spouse:
- Guarantor:
- Legal Guardian:
- Other Contact:
- Email Home: mikemcnair@msn.com
- Email Work:
- Insurance Note: Delta Dental of VA
- Web ID: MMcNair
- Web Password: *

Prefer Correspondence via:

- ☐ Mail ☐ Fax ☐ Email (H) ☐ Email (W)
- ☐ Patient will call to set up appointment
- ☐ Need to call patient to set up appointment
- ☐ Do Not Schedule
- ☐ Deceased
- ☐ Lecture
- ☐ Board
- ☐ Export
- ☐ Implant

Remove Patient

Control Table Information

Practice Information Application Settings Network Settings Workstation Settings

Digital Radiography
Radiography Software: Trophy/Kodak
Radiography Path/Source: C:/Program Files/BreezeSys\DSLR Remote Pro

Scan Twain Settings
☒ Default ☐ Advanced ☐ Remote

Remote Camera
Software: Breeze Systems
Remote Camera Software Path: c:\Program Files\BreezeSys\DSLR Remote Pro
☐ Autosave Camera Image ☐ Autoload Remote Software
☐ Auto Archive ☒ Include Subfolders

Inactivity Duration Before Logoff
Minutes: 20

Credit Card Settings
Credit Card Company: Worldpay Connection Method: Com Port: Pin Pad Type: ☐ Canadian EMV

Image Location
Flash Card Reader Path: C:/Program Files/Brown & Herbranson\3D Tooth Atlas

Patient Privacy
☒ Open Application Showing Calendar First

Tooth Models
Tooth Model Path: C:/Program Files/Brown & Herbranson\3D Tooth Atlas

Digital Signature
Signature Type: Leadtools (Default)

Notes Encryption
☒ Allow Notes Encryption

Default Image Category
Leadtools (Default)

Windows Placement
Reset ☒ Allow TDO to be Full Screen

Close

Record: 14 of 3 No Filter Search

- Under Digital Radiography, select Trophy/Kodak. Change the path to C:/TW or C:/program files/TW, depending on where the TW is located.
- C:\Program Files (x86)\Kodak\KODAK Dental Imaging Software -is an additional path that should be added as a possible path.
- Then, when you select the tooth icon, it will open Carestream Imaging.
- You can use that software to acquire an image.
- Once you have your image in Carestream image, copy the image, then paste the image into the acquire window and save it in TDO normally.

Online URL: <https://kb.tdo4endo.com/article.php?id=860>